



Nicolet® Clinical EEG Modular Neurodiagnostic System

Technical Specifications

General Specifications

Isolated Power Supply 110/115 or 230 VAC \pm 10% input, 50 - 60 Hz 595 VA primary; 500 VA secondary Output voltage = input voltage

Dimensions Weight................ Unibody cart approx. 68kg (150 lbs.) (depending on model of printer) Operating Environment (in use) Non-Operating Environment (in storage)

Computer

The Nicolet Clinical EEG system operates with a desktop computer. Please contact your Natus representative for the latest computer specifications.

HP DeskJet printer (black, white, and color)

Network

10/100/1000 Mb Ethernet (standard) HL7 compatible through NicVue Connect interface

EEG software (with v32/v44 Amplifier)

EEG Display Sensitivity...... 10, 20, 30, 50, 70, 100, 150, 200, 300, 500, 700, 1000, 2000, 5000 μ V/cm 1, 2, 3, 5, 7, 10, 15, 20, 30, 50, 70, 100, 200, 500 μV/mm Off, 10, 15, 25, 30, 35, 40, 50, 60, 70, 100, 150, 200, 300, 500, 1000, 1500 Hz Low Filters Off, 0.053, 0.16, 0.3, 0.5, 1, 1.6, 2, 3, 5 Hz 0.2, 0.33, 0.5, 0.625, 1, 2, 3.3, 6.2, 18.9 seconds

Notch FilterOff, 50/60 Hz



v32 Amplifier

Analog/Digital Converter	16 bits
ADC Resolution Voltage	= 0.153 μV
DC Offset Tolerance	± 340 mV
Channels (Inputs) 32 EEG, configurable as bipolar AC (2	24-32), 1 configurable as DC (32)
Maximum Input Range	\pm 5 mV
Bandwidth	0.053 - 500 Hz
Noise	< 1.5µV pk-pk @ 0.1 - 100 Hz
Input Impedance	> 100 M Ω (common mode)
CMRR at Patient Inputs > 115 dB @ $50 - 60$ Hz, with	active patient ground connected
Channel Crosstalk	
Amplifier Sample Rate (under software control)	125, 250, 500, 1000, 2000
Calibration Square wave, 1, 5, 10, 20 sec period,	10, 50, 100, 1000 μV amplitude
Input Bias Current	< 5 nA
Anti-Aliasing Filter Cut Off Frequency	500 Hz
Differential Input Impedance	40 ΜΩ
Interface to Amplifier	Ethernet
Built-in Impedance and Display	
Headbox	. Optional; no impedance display
Additional Ports	
 Isolated SpO₂ with X-Pod 	
Photic output	
 Isolated patient event button 	
Channel Hardware Gain	410
Deblock	Yes
Auxiliary Inputs	
1 Hi-level, non-isolated input for connection of external devices (e.g. CO2 monitors, etc.)	
Analog/Digital Converter	16 bits
Maximum Input Range	
ADC Resolution	76.3 μV

www.natus.com

Nicolet Clinical EEG

v44 Amplifier

System Configurations

Sleep, EEG, ICU monitoring and LTM

OR and non-OR applications

Cart mount and wall mount options
Analog/Digital Converter
ADC Resolution Voltage = $0.153 \mu V$
DC Offset Tolerance \pm 900 mV
Channels (AC Inputs)32 EEG, (9 configurable as bipolar (24-32) AC)
12 non-isolated DC inputs (\pm 5 V, BW = 100 Hz)
Maximum Input Range ± 5 mV
Bandwidth
Noise< 1.5μV p-p @ 0.1 - 100 Hz (except channels 31, 32 and
OR channels < 2uV p-p @ 0.1 - 100 Hz)
Input Impedance > 100 M Ω (common mode)
CMRR at Patient Inputs > 115 dB @ $50 - 60$ Hz, with active patient ground connected
(except channels 31, 32 and OR channels > 100 dB @ 50-60 Hz with RLD)
Channel Crosstalk< -40 dB
Amplifier Sample Rate (under software control) 125, 250, 500, 1000, 2000
Calibration Square wave, 1, 5, 10, 20 sec period, 10, 50, 100, 1000 μV amplitude
nput Bias Current< 5 nA
Anti-Aliasing Filter Cut Off Frequency500 Hz
Differential Input Impedance
nterface to Amplifier Ethernet
Channel Hardware Gain

Deblock Yes

Integrated SpO₂

Channels (DC Inputs): 12 non-isolated

- Analog/Digital converter: 16 bits
- Maximum input range: ± 5V
- ADC resolution: 153 μV
- Bandwidth: DC 100 Hz

Additional Ports

- Panasonic camera control port on amplifier
- Isolated SpO₂
- Isolated patient event button
- Photic output

Headboxes

v44 requires one of the following:

- Clinical headbox with built in impedance and display
- Clinical headbox with head cap adapter and built in impedance and display
- OR headbox

Interfaces with the Nicolet Ambulatory Module. See Nicolet Ambulatory EEG Specifications Sheet #169-435300.

Quality Standards

 \bigoplus

Manufactured, designed, developed and marketed under ISO 13485 certified quality system

Compliance/Regulatory Standards

Designed, tested, manufactured and certified to meet the following domestic (USA),

Canadian, European and International Standards:

UL 60601-1 Medical Electrical Safety Standard (USA)

CAN/CSA-C22.2 no. 601.1-M90 Medical Electrical Safety Standard (Canada)

EN/IEC 60601-1 Medical Electrical Safety of Medical Equipment (International and Europe)

IEC 60601-2-26 Particular Safety of electroencephalographs equipment

IEC 60601-2-40 Particular Safety of electromyography and evoked response equipment

EN 60601-1-2 Collateral safety standard for EMC

European Community (CE Mark)

Medical Device Directive (MDD) product certified to comply to EC Directive 93/42/EEC Patient IsolationBF



Natus Neurology Incorporated

3150 Pleasant View Road Middleton, WI 53562 USA 1-800-356-0007

1-608-829-8500 Fax: 1-608-829-8709 www.natus.com

Specifications subject to change without notice.

©2016 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates.

P/N 169-432801 Rev10



